

# MAGNUM 3500 PRO

## PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

THE ONLY GLOBUS MAGNETOTHERAPY DEVICE WITH FOUR CHANNELS THAT ALLOWS TO TREAT DIFFERENT BODY PARTS OR SUBJECTS AT THE SAME TIME. MAGNUM 3500 PRO IS EQUIPPED WITH A WIDE RANGE OF ACCESSORIES FOR ALL TREATMENT NEEDS.

### TECHNICAL SPECIFICATIONS

#### Display

Size of the visible area 2.6"

#### Channels

4 independent

#### Intensity can be emitted with flexible diffuser

- 200 Gauss on average for channel
- 800 Gauss on average on the whole (for 4 channels)

#### Frequency

5-200Hz

#### Duty Cycle

5-50%

#### Free programming

Possibility of adding up to 10 custom programs

#### 2+2 function

Possibility of running 2 different programs at the same time, using 2 channels each

#### Power supply


- Mains
- Rechargeable batteries (not included)

#### Size | Weight

170x220x60 mm | 1152 gr

### STANDARD EQUIPMENT

- 1 Suitcase
- 1 MAGNUM 3500 PRO magnetotherapy device
- 2 Flexible solenoids 330x110 mm
- 2 Soft solenoids 120x125 mm
- 1 Soft mat 4 250x250 mm
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating manual



CE  
0476



- EASY PROGRAM
- 2+2 MODE
- LAST 10
- PROGRAMMABLE
- RENTAL OPTION
- RUN TIME / EDIT PARAMETERS
- MAINS SUPPLY

**GLOBUS**  
ITALIAN EXCELLENCE

# MAGNUM 3500 PRO

## PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

### PROGRAM LIST

Rev. En\_01\_2020

#### ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthritis

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

#### PAIN

Algodystrophy

Back pain

Joint pain

#### INFLAMMATIONS

Anti-inflammatory

Tendinitis

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

#### ORTHOPEDICS

Post ACL surgery

Cartilage lesions

Knee prosthesis

Hip prosthesis

Fusion

Osteotomy

#### BONE PATOLOGIES

Fractures

Scaphoid fracture

Wrist fracture

Shoulder fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteobecrisis

Osteoporosis

Pseudoarthrosis

Bone edemas

#### TRAUMAS/EDEMAS/WOUNDS

Ankle sprain

Edemas - Contusion

Ulcers

#### NIGHT PROGRAMS

Night program 15Hz

Night program 50Hz

Night program 75Hz

#### DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multi-frequency program

Program 5Hz

Program 10Hz

Program 15Hz

Program 20Hz

Program 25Hz

Program 30Hz

Program 35Hz

Program 40Hz

Program 45Hz

Program 50Hz

Program 55Hz

Program 60Hz

Program 65Hz

Program 70Hz

Program 75Hz

Program 80Hz

Program 85Hz

Program 90Hz

Program 95Hz

Program 100Hz

Program 120Hz

Program 140Hz

Program 160Hz

Program 180Hz

Program 200Hz